
Helpful Email References: SCANS By DATE

Hi Diana,

We have a question is it possible when you want to check the KOD form several date's that the computer say the intestine when you look that up that the first date is on top and the last scan is at the end with everything in between?

No you see that say first the 1th of December than 28th of November than 30th of November it's not correct on date. Or is this possible to do with sorting on the name and the data at the same time?

See the enclosed picture

Many thanks and all the best,

29.11.18	12:14:21		BLOOD CELLS	S		1,548	43.43
1.12.18	12:49:10		BLOOD CELLS	S		1,539	64.57
2.12.18	13:32:28		BLOOD CELLS	S		1,545	58.28
28.11.18	13:43:25		BLOOD CELLS	S		1,544	54.86
30.11.18	12:28:47		BLOOD CELLS	S		1,516	62.86

You can already do this.

1. type in the name of the organ you are looking for (entire thing)
2. You will have them there, click on "DATE" and everything will get in order

DATE	TIME	THUMB	NAME	PF	Sys	COLOR	KOD
11/1/2018	08:58:13		STOMACH WALL	S		Nor	1,351
11/6/2018	11:22:06		STOMACH WALL	S		Acu	1,287
11/26/2018	16:50:27		STOMACH WALL	S		Acu	1,301
11/27/2018	11:43:18		STOMACH WALL	S		Acu	1,301
11/30/2018	08:43:28		STOMACH WALL	S		Acu	1,191